


For chronic and acute wounds

LAVANID®





SERAG-WIESSNER specialises in moist wound management and is your expert in this field, thanks to a perfectly matched product range with a wide selection of various product containers and dosage forms.

The LAVANID® product range is suited for cleansing and moistening both chronic and acute wounds as well as for intraoperative use.

wound irrigation solution based on Ringer's solution with 0.04% polyhexanide as preservative

good tolerability

Fields of application of LAVANID® Wound Irrigation Solutions:

acute and chronic wounds
intraoperative wounds
burns
NPWT (negative pressure wound therapy)





LAVANID® - Wound Irrigation Solution

LAVANID® is a sterile, isotonic wound irrigation solution based on Ringer's solution with 0.04% polyhexanide as preservative. Polyhexanide prevents

the growth of microorganisms in wound dressings and in the solution.

The isotonic solution prevents swelling of the skin and reduces an irritation.

Wound management products preserved with polyhexanide, such as LAVANID®, are the treatment of choice for poorly healing chronic wounds and sensitive wounds, including e.g. burns.

painless application

reduces inflammation signs

prevents biofilm formation



Pouring out the solution with screw top.



Bottle with inserted spike ready for aseptic removal.



There are no restrictions on duration or frequency of use. A post-treatment or rinsing off is not necessary.

Method of application

LAVANID® - Wound Irrigation Solution:
The stopper allows aseptic removal of the irrigation solution, e.g. with a spike. Alternatively, the screw cap can be removed to pour out the solution.

Composition

Ringer's solution (isotonic electrolyte solution consisting of sodium chloride, potassium chloride, calcium chloride, and water for injection purposes), polyhexanide 0.04%.

Shelf life:

2 months after initial opening

Delivery programme - LAVANID® Wound Irrigation Solution

Pack sizes	PZN	REF
20 x 125 ml	18660576	018024
10 x 250 ml	18660582	018025
10 x 500 ml	18771028	018026
6 x 1000 ml	18660599	018027



LAVANID® - Wound Gel

LAVANID® - Wound Gel is a sterile hydrogel based on Ringer's solution and 0.04% polyhexanide as preservative. Polyhexanide prevents the

good flow properties

prevents biofilm formation

good tolerability

growth of microorganisms in wound dressings and in the gel. The hydrogel exhibits good flow properties, for example when applied in wound cavities or pockets.

Application of LAVANID® - Wound Gel/ LAVANID® - Wound Gel V+

First, clean the wound with LAVANID® Wound Irrigation Solution. Apply the wound gel directly to the wound or wound dressing, or use a sterile spatula, then cover the wound with a dressing. The wound gel remains in place on the wound until the next dressing change. Any residual wound gel can be rinsed off with LAVANID® Wound Irrigation Solution.



Fields of application of LAVANID® - Wound Gel/LAVANID® - Wound Gel V+

moistening and cleansing acute and chronic wounds

burns

loosening wound coatings (autolytic debridement)

prevents biofilm formation



LAVANID® - Wound Gel V+

LAVANID® - Wound Gel V+ is a sterile hydrogel based on Ringer's solution and 0.04% polyhexanide as preservative. Polyhexanide prevents the growth of microorganisms in wound dressings and in the gel.

Product-related pain is not expected during application. Both hydrogels have a high water content (approximately 90%) and support autolytic debridement.

A clinical study on soft tissue wounds showed that the wound gel was well tolerated, and that no local irritation or pain is to be expected.

Composition

Ringer's solution (consisting of sodium chloride, potassium chloride, calcium chloride, and water for injection purposes), glycerin, hydroxyethyl cellulose, polyhexanide 0.04%.

better adhesion compared with conventional wound gel

prevents biofilm formation

good tolerability

Delivery programme - LAVANID® Wound Gel

Pack sizes	PZN	REF
1 x 15 g tube	17817136	015275
1 x 40 g tube	13946724	015273
1 x 100 g tube	08837192	015274
2 x 10 g syringe	06122661	015272
12 x 10 g syringe	02078785	015271

Delivery programme - LAVANID® Wound Gel V+

Pack sizes	PZN	REF
1 x 40 g tube	10813543	015773
1 x 100 g tube	10329497	015774
2 x 10 g syringe	13576919	015772
12 x 10 g syringe	13576925	015771

Anticipated risks associated with the intended use of the products are specified in the instructions for use. These are available on the SERAG-WIESSNER website at www.serag-wiessner.de.



**SERAG
WIESSNER**

SERAG-WIESSNER GmbH & Co. KG
Zum Kugelfang 8 - 12
95119 Naila/Germany



+ 49 9282 937-0



+ 49 9282 937-9369

Export Department:



+ 49 9282 937-230



+ 49 9282 937-9785



info@serag-wiessner.de



www.serag-wiessner.de

